

Medimate MiniLab







> 12 years research & development



2000: Initiative

Psychiatrist Pieternel Kölling Prof. Albert van den Berg

2001: PhD Research

2006: Founding Medimate

2007: Advisory Board





Result

Lithium measurements available any place, any time









Finger stick

Droplet on lab-chip

Lab-chip in device

Li result



Result

Everyone can now measure Lithium:

- several minutes
- everywhere
- only one droplet blood (capillary, serum and venous)
- plug and play
- accurate





To improve quality of live

Traditional

several minutes

MEDIMATE

> 1 hour Hospital

Home

KISS

Complex

Patient

Highly educated

Tube of blood (Serum)

1 droplet of blood

No

Yes

Patient Empowerment:

Analyze:

Where:

User:

Sample:

Hardware:





Selfmonitoring

- pregnancy
- difficult to go to laboratory
- difficult to match opening hours lab
- travelling
- changes in comedication
- instable





Docters desk

- toxication?
- mentally instable
- starting lithium
- elderly
- pregnancy
- fear
- therapy ontrouw





Specifications

Reach:

0,2 - 7,0 mmol

Outside: < 0,2 of > 7,0

Expiry date:

9 months





Tolerance

Allowable error

CLIA, WSLH, AAB: +/- 20% or 0.3mmol/L

NYS: +/- 15% or 0.3mmol/L

CAP: 20% or 0.3mEq/L

RCPA: 0.2mmol/L


CFX: 10%


+/- 0.1% or 0.1 mmol/L






Accuracy












Fitted for data transfer







Cost

- Lithium lab-chip € 18,50 excl. VAT
- Multireader approx. € 1.250 excl. VAT





Rollout

2012 Validation:

MST: clinical - chemistry

Dimence: accuracy and functionality users

2013 Trimbos Intitute: in use by user

75 patients, 12 month, Lithium MiniLab at home

cost / tehrapy effectivity / well being

2013 **GGZ** Friesland: in use within organisation





Rollout







Rollout

The Netherlands (with Health Insurance Companies)
Spain (with Francesc Colom)
???

2013 Q4 - FDA for USA

2014 USA ???





Help

Introduction in your country:

- Opnion leaders
- Psychiatry
- Pharmacology
- Clinical chemistry
- National advocacy networks





+31(0)88 66 55 222 of +31(0)6 1654 2848
huub.maas@medimate.com
marian.tellegen@medimate.com



